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(GB)

- For easy blackhead removal
- Helps to unclog & reduce pores
- Rejuvenates and refreshes skin

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(GB) The multifunctional Pore cleansing System is designed to help remove blackheads and increase blood circulation, while ridding skin of dirt and impurities. Features an ergonomic handle for a comfortable and easy use.

An innovative, battery-operated device to help you cleanse and clear your skin!

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**(GB) HOW TO USE:**

- Prior to first use, insert 2 x AA batteries following the direction for(+/-)
- Test the device on your arm to determine how skin will react to the treatment
- Wash skin with lukewarm water to open up pores before using the device to ensure easier removal
- Use only on T-Zone areas(forehead, nose and chin)
- Gently glide the device over face using the appropriate nozzle
- Do not repeat or treat one area for too long because bruising/irritation may occur
- Follow treatment using desired toner, serum or moisturizer

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**(GB) REPLACING THE FILTER:**

1. Remove the nozzle from the device
2. Remove filter cover and disinfect with alcohol
3. Remove filter with tweezers and/or small pliers to disinfect the area
4. Replace with a new filter and close the filter cover

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**DIAMOND HEAD:** Use to gently exfoliate away dead skin cells, make-up and impurities to reveal fresher, smoother skin

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**LARGE CIRCULAR HEAD:** Use to effectively remove blackheads, increase blood circulation and rejuvenate skin

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**SMALL CIRCULAR HEAD:** Use on more precise and sensitive areas such as under the eyes

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**OVAL HEAD:** Use to firm skin and reduce the appearance of fine lines and wrinkles

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**(GB) CARE INSTRUCTIONS:** After each use, remove the nozzle from the device and disinfect with alcohol. Lay flat to air dry. Ensure that you replace filter sponge after 4-6 uses.

**CAUTION:** Do not use in one area for too long because bruising may occur. Do not repeat areas already treated. Use only as directed. Only use on T-zone areas. Never use on sensitive, irritated, blemished or sunburned skin. Discontinue use if skin becomes red, swollen or itchy and consult a doctor if symptoms persist after use. For external use only. Avoid contact with eyes. Keep out of reach of children.

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#### **GB- CAUTION BATTERY:**

1. Non-rechargeable batteries are not to be recharged
2. Never mix different types of batteries
3. Empty batteries are to be removed from the item. Do not mix them with the new ones.
4. Never short-circuit the contacts. Follow the battery polarity indications.
5. Change the batteries if they do not work properly or they have leak.
6. Never put the batteries near the fire or heat sources, nor throw them into the fire.
7. Please wash immediately if the skin or clothes touched the liquid from the battery.
8. Changing of batteries to be done by adults. Keep away from children.
9. In case of swallowing seek medical attention immediately!